# **TRANSFECTION**

WORLD BEST PERFORMANCES

**WELGENE**

---

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# WelFect-EX™ Transfection Reagent

## WelFect-EX™ PLUS Transfection Reagent

### WelFect-EX™ GOLD Transfection Reagent

For Broad Cell Lines, Up to 100% Transfection Efficiency

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

Company A

#### ● Transfection Efficiency 2

- + Cell Line: HeLa
- + Photo: Fluorescent microscopy (X 100)  
(DNA = 0.4 ug/well, pcDNA – GFP)

WelFect-EX™ GOLD

WelFect-EX™ PLUS

WelFect-EX™

Company A

#### ● Ordering Information

Product Name	Cat. No.	Size
WelFect-EX™ PLUS Transfection Reagent	TR 003-01	0.5 ml (300 Transfections)
WelFect-EX™ PLUS Transfection Reagent	TR 003-02	1 ml (600 Transfections)
WelFect-EX™ PLUS Transfection Reagent	TR 003-03	1 ml X 5ea (3,000 Transfections)
WelFect-EX™ GOLD Transfection Reagent	TR 008-01	0.5 ml (300 Transfections)
WelFect-EX™ GOLD Transfection Reagent	TR 008-02	1 ml (600 Transfections)
WelFect-EX™ GOLD Transfection Reagent	TR 008-03	1 ml X 5ea (3,000 Transfections)
WelFect-EX™ Transfection Reagent	TR 006-01	0.5 ml (300 Transfections)
WelFect-EX™ Transfection Reagent	TR 006-02	1 ml (600 Transfections)
WelFect-EX™ Transfection Reagent	TR 006-03	1 ml X 5ea (3,000 Transfections)

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# WelFect-Q™ GOLD Transfection Reagent

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





## ● Transfection Efficiency

## ● Ordering Information

Product Name	Cat. No.	Size
WelFect-Q™ PLUS Transfection Reagent	TR 002-01	0.5 ml (300 Transfections)
WelFect-Q™ PLUS Transfection Reagent	TR 002-02	1 ml (600 Transfections)
WelFect-Q™ PLUS Transfection Reagent	TR 002-03	1 ml X 5ea (3,000 Transfections)
WelFect-Q™ GOLD Transfection Reagent	TR 009-01	0.5 ml (300 Transfections)
WelFect-Q™ GOLD Transfection Reagent	TR 009-02	1 ml (600 Transfections)
WelFect-Q™ GOLD Transfection Reagent	TR 009-03	1 ml X 5ea (3,000 Transfections)

# ENHANCER-Q™ Transfection Enhancing Reagent


# ENHANCER-GOLD™ Transfection Enhancing Reagent

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## Enhancing of Transfection Efficiency

HeLa  
pcDNA 3.0 – GFP  
Fluorescent  
microscopy

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
## Oncogene 30, 3345–3359

### Transfection

Colon cancer cells ( $5 \times 10^4$  cells per well) were plated in 24-well plates and transiently transfected with 0.4  $\mu\text{g}$  of the empty vector or the constitutively activated full-length IL-32 $\gamma$  plasmid per well, using a mixture of plasmid and the **WelFect-EX PLUS** reagent in OPTI-MEN, according to the manufacturer's specification (WelGENE, Seoul, Korea).

### Cell viability

To determine viable cell numbers, the colon cancer cells were seeded onto 24-well plates ( $5 \times 10^4$  cells per well). The cells were trypsinized, pelleted by centrifugation for 5 min at 1500 r.p.m., resuspended in 10 ml of PBS and 0.1ml of 0.2% Trypan blue was added to the tumor cell suspension in each solution (0.9 ml each). Subsequently, a drop of suspension was placed in a Neubauer chamber, and the living cancer cells were counted. Cells that showed signs of Trypan blue uptake were considered to be dead, whereas those that excluded Trypan blue were considered to be viable. Each assay was carried out in triplicate.

## J. Biol. Chem. 2010 285: 36410–36419.

### Gene Regulation

Ji Hae Jun, Won-Joon Yoon, Sang-Beom Seo, Kyung-Mi Woo, Gwan-Shik Kim, Hyun-Mo Ryoo, and Jeong-Hwa Baek

#### BMP2-activated Erk/MAP Kinase Stabilizes Runx2 by Increasing p300 Levels and Histone Acetyltransferase Activity

J. Biol. Chem. 2010 285: 36410-36419. First Published on September 17, 2010, doi:10.1074/jbc.M110.142307

...Lipofectamine 2000 reagent was purchased from Invitrogen and **WelFect-Ex™ PLUS** was from WelGENE (Seoul, Korea). The luciferase assay...transiently transfected with pcDNA or MEK1 expression vectors using WelFECT. After overnight recovery, the cells were incubated in the...

## J. Biol. Chem. 2011 286: 7439–7456.

### Signal Transduction

Hyun-A Seong, Haiyoung Jung, and Hyunjung Ha

#### Murine Protein Serine/Threonine Kinase 38 Stimulates TGF- $\beta$ Signaling in a Kinase-dependent Manner via Direct Phosphorylation of Smad Proteins

J. Biol. Chem. 2010 285: 30959-30970. First Published on July 21, 2010, doi:10.1074/jbc.M110.138370

...were transfected with the p3TP-Lux or p21-Luc reporter plasmids, along with each expression vector as indicated, using **WelFect-Ex™ PLUS**. The cells were lysed, and luciferase activity was detected using the dual luciferase assay kit (Promega, Madison...

## FEBS Letters 584 (2010) 3923–3928

### 2. Materials and methods

#### 2.1. Plasmids

Full-length mouse cDNA of PTPRO was ligated with pCS2-HA or pCMV5 vector to make a PTPRO-HA or PTPRO-Flag construct, respectively. To make PTPRO(ECD)-Flag or PTPRO( $\Delta$ Ex)-Flag, PCR products containing the region from the 5' to 3' end of the extracellular domain or from the transmembrane domain to the 3' end of the mouse PTPRO cDNA were cloned into pCMV5 vector, respectively. For PTPRO-ECD(1-3)-Fc, a PCR product containing the region from the 5' end to the third FN III repeats of the extracellular domain was cloned into pRK-Fc vector. To make Ptpo-Luc, about 2 kb from the transcription initiation site of *ptpo* was amplified by PCR and ligated into pGL3-basic vector.

#### 2.2. Cell culture and transient transfection

HEK293T, L929 and CaCO<sub>2</sub> cells were maintained in Dulbecco's Modified Eagle's Medium (DMEM, Gibco BRL) supplemented with 10% Fetal Bovine Serum (FBS) and antibiotics. For transient transfection experiments, each plasmid DNA was transfected by the calcium phosphate precipitation method. For L929 cells, **Welfect-Ex™** was used according to the manufacturer's instruction (Welgene). The amount of DNA in each transfection was held constant by the addition of an appropriate amount of empty expression vector.

## Cellular Signalling 21 (2009) 43 – 51

### 2. Materials and methods

#### 2.1. Cell cultures and transient transfection

C57mg, HEK 293T and MCF7 cells were grown in Dulbecco's modified Eagle medium (DMEM, Life Technologies) supplemented with 10% fetal bovine serum (FBS, Cambrex Bio Science Walkersville, Inc.) and 1X antibiotic-antimycotic (catalog number 15240; Invitrogen) in humidified 6% CO<sub>2</sub>. Cells were transiently transfected using either the calcium phosphate precipitation method or **WelFect-Ex™ PLUS** transfection reagent (WelGENE), in accordance with the manufacturer's instructions. Leptomycin B (LMB) was purchased from Sigma and cells were treated with 50 ng/ml for 8 h before immunofluorescence analysis.

#### 2.2. Reverse transcriptase-PCR

Total RNAs isolated from cell lines, different tissues and embryos of CD 1(ICR) mice were reverse transcribed using random primers

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